



Soft Newborn Z-Band Wristband Solutions

Gentle enough for a newborn's delicate skin

Wristbanding newborns is a sensitive issue. While hospitals go to great lengths to properly identify newborns, clinicians are also concerned with the comfort and durability of the bands. Zebra kept these challenges in mind when creating a complete line of wristbands for newborns.

Barcoded wristbands are the foundation for positive patient identification and the prevention of serious medical errors. Zebra thermal printers produce crisp, clear print quality on all wristbands, so barcodes scan quickly and reliably. Zebra offers several soft, comfortable newborn wristbands to meet most hospital's needs.


Zebra's ZD510-HC Patient I.D. Solution


The reliable ZD510-HC direct thermal printer combined with easy-to-load cartridges of Z-Band wristbands in adult, pediatric or newborn sizes allow patient I.D. protocols to be easily implemented. Pop the wristband cartridge into the ZD510-HC printer and print durable, individual wristbands on demand. The ZD510-HC printer detects the wristband size and automatically calibrates its settings for optimal print quality.



FACT SHEET




SOFT NEWBORN Z-BAND WRISTBAND SOLUTIONS

	Z-Band Newborn Wristbands	Z-Band Newborn Soft Wristbands
	<p>Flexible polypropylene material with a soft-feel coating</p> <ul style="list-style-type: none"> Part Number: 10035439K 	<p>Soft nylon backing for additional comfort</p> <ul style="list-style-type: none"> Part Number: 10031289K
	<ul style="list-style-type: none"> Adjustable band fits the smallest patients print area Flat print area ensures ease of scanning 0.75 in. x 7.69 in./19mm x 195mm band - 0.75 in. x 2 in./19 mm x 51 mm print area Patented antimicrobial coating protects the wristband Pressure sensitive, adhesive tab closure for quick and easy securement 	

Z-Band Fusion ID Tags	
	<p>Zebra's Z-Band Fusion ID Tags feature a self-laminating design that protects the image, making them perfect for use in the NICU where the average newborn stay exceeds 14 days.</p> <ul style="list-style-type: none"> Can be used in conjunction with the ComfyCuff® foam wristband, SwaddleBand and Vinyl Newborn Wristbands 1.75 in. x 0.875 in./44.5 mm x 22.23 mm (1 in./25 mm core) Latex-free Part Number: FID-BABY-T-1-500

Newborn Wristband Accessories

Zebra's ComfyCuff family of wristband accessories when paired with Z-Band Newborn wristbands, offers clinicians a way to identify newborns without having to unwrap them.

<p>ComfyCuff</p>  <p>The soft, foam cuff features an adjustable Velcro closure ensuring newborn comfort. It can be used in conjunction with Z-Band Fusion ID tags.</p> <ul style="list-style-type: none"> 5 in. x 1.25 in./27 mm x 32 mm Latex-free Part Number: LB2-BABY-P 	<p>ComfyCuff Security Seal</p>  <p>Label securely sticks to the ComfyCuff and tears away to indicate tampering.</p> <ul style="list-style-type: none"> 2.875 in. x 0.4375 in./73 mm x 11 mm Latex-free Part Number: LB2-BABY-CCS
<p>ComfyCuff Newborn Swaddle Band</p>  <p>Tyvek® band creates an extension between the ComfyCuff and Z-Band Direct Soft Newborn Wristband or Z-Band Fusion ID Tag to make identification of swaddled babies quick and easy.</p> <ul style="list-style-type: none"> 9 in./229 mm long Latex Free Part Number: LB2-BABY-SWB9 	<p>Vinyl Newborn Wristbands</p>  <p>These latex-free, soft, flexible vinyl wristbands are for use with Z-Band Fusion newborn ID tags. The narrow width allows some models of "Abduction Prevention Systems" to be attached. They are available in clip or self-locking closure. The self-locking closure wristband enables size adjustment with a gentle tug as the newborn loses mass around the wrists and ankles within the first 24 hours of birth.</p> <ul style="list-style-type: none"> Clip-Closure Part Number: PS-BABY Self-Locking Closure Part Number: PS-NEWBORN-ADJ

Specifications subject to change without notice.

For more information about Z-Band Newborn Wristbands, visit www.zebra.com/wristbands



NA and Corporate Headquarters
+1 800 423 0442
inquiry4@zebra.com

Asia-Pacific Headquarters
+65 6858 0722
contact.apac@zebra.com

EMEA Headquarters
zebra.com/locations
contact.emea@zebra.com

Latin America Headquarters
+1 786 245 3934
la.contactme@zebra.com